

## **Viniferamine® Wound Hydrogel: The Occupied Wound Interface and a New Approach to Biofilm Prevention**

Darlene E. McCord, PhD

### **Introduction: The Wound Surface as a Physical Environment**

One of the persistent challenges in wound care is biofilm. Much of the discussion surrounding biofilm has understandably focused on how microorganisms are killed after they enter the wound. Anti-microbials, antiseptics, debridement, cleansing, and other interventions are commonly directed toward reducing an established microbial burden.

Viniferamine® Wound Hydrogel suggests an additional question: What if the physical architecture of the wound interface could make it more difficult for biofilm to establish itself in the first place?

This concept begins with the hydrogel's unusual physical behavior. When applied at approximately 60–65°F (15.5–18.3°C), the material has relatively low viscosity. It can therefore flow across an irregular wound bed and enter its contours, depressions, and surface variations. As the hydrogel quickly warms upon contact with the wound, its viscosity increases, forming an intimate, conformal mold of the wound surface.

The importance of this behavior may extend considerably beyond ease of application. The resulting structure can be viewed as an occupied wound interface—a wound surface in which the hydrogel intimately conforms to the available wound topography rather than simply resting above it. That distinction could have important implications for biofilm.

### **Biofilm Requires More Than Microorganisms**

Biofilm formation is not simply the presence of bacteria in a wound. It is an organized biological process. Microorganisms must first reach the wound surface. They must encounter an environment in which attachment is possible, adhere to tissue or another surface, begin multiplying, and produce extracellular polymeric substances that contribute to the development of the biofilm matrix. This means that the earliest stages of biofilm development depend upon an interface. The physical character of that interface may therefore matter.

An irregular wound covered by a dressing does not necessarily mean that every portion of the wound is actually in intimate contact with the dressing. Wound topography can create valleys, crevices, undermined areas, and other microscopic or macroscopic spaces beneath or within the dressing interface.

Exudate may accumulate within these spaces. Microorganisms carried by wound fluid may enter the wound. If organisms encounter exposed wound surfaces within a protected microenvironment, conditions may be conducive to attachment and subsequent colonization. This leads to an important hypothesis: Gaps in the wound interface may represent a biological opportunity.

## **The Concept of an Occupied Wound Interface**

Viniferamine Wound Hydrogel is designed to behave differently. Its initial low viscosity allows the material to conform to the wound architecture rather than forcing the wound to conform to the geometry of a conventional dressing.

Once intimate contact has been established, warming increases viscosity and helps preserve that conformity. The result is an intimate mold of the wound surface. Instead of placing a relatively planar material over an irregular biological surface, the hydrogel becomes geometrically matched to that surface. This may fundamentally change the wound-dressing interface. The goal is not merely to cover the wound. The goal is to occupy it.

### **No Gaps. No Voids. Fewer Places to Anchor.**

This distinction may provide a new way of thinking about biofilm prevention.

Microorganisms cannot establish a mature biofilm without first interacting with a surface. Attachment is one of the earliest critical events in the progression toward biofilm formation. A conformal hydrogel that intimately occupies the wound surface may reduce the unfilled spaces available for pooled exudate and may limit exposed tissue interfaces within those spaces.

Conceptually, the sequence becomes: Flow → Conform → Occupy → Maintain Contact → Reduce Opportunities for Attachment.

This is fundamentally different from relying exclusively upon antimicrobial activity after microorganisms have attached. It introduces a physical component to biofilm management. The hydrogel does not merely deliver something to the wound. Its physical presence at the wound interface may be important in itself.

## **A Bioidentical Mold**

The term bioidentical mold provides a useful description of this phenomenon. Every wound has a unique three-dimensional architecture. Even wounds that appear relatively smooth macroscopically contain irregularities at smaller scales.

A sufficiently conformable hydrogel can adapt to this architecture during application. Rather than bridging from one elevated region of tissue to another and leaving a void beneath, the material can flow downward and establish closer contact with the underlying wound surface.

As it subsequently increases in viscosity, that geometry is retained. The wound, therefore, effectively determines the shape of its dressing interface.

The treatment adapts to the wound rather than requiring the wound to adapt to the dressing. For wounds with complex or irregular surfaces, this could be especially valuable.

## **The Occupied Interface Hypothesis**

This concept can be expressed as a testable scientific hypothesis:

A highly conformal wound hydrogel that forms an intimate mold of the wound surface will reduce interfacial voids and exposed microenvironments, thereby reducing opportunities for microbial attachment and subsequent biofilm establishment compared with less-conformal wound interfaces.

This hypothesis is important because it does not rely on the claim that physical conformity alone kills microorganisms. Instead, it proposes that physical occupation of the wound interface may interfere with one of the prerequisites for biofilm formation: access to an attachment environment.

Antimicrobial activity and physical interface occupation could therefore represent

complementary mechanisms. One addresses the microorganism. The other addresses the environment in which the microorganism can establish itself.

### **From Biofilm Treatment to Biofilm Prevention**

This distinction changes the way we think about wound biofilm. Traditional thinking often begins after colonization has occurred:

- Biofilm forms → Biofilm matures → Treatment attempts to disrupt or kill it.
- The occupied-interface concept moves attention toward an earlier point:
- Wound is treated → Interface is occupied → Attachment opportunities are reduced → Biofilm establishment may become more difficult.

Instead of asking only, "How effectively can we destroy biofilm?" we can also ask, "How can we design the wound environment so that biofilm has fewer opportunities to establish itself?"

### **The Importance of Exudate Management**

This concept does not eliminate the importance of exudate management. A conformal hydrogel should not simply trap uncontrolled fluid against the wound. The wound-care system must provide an appropriate pathway for excess exudate to move away from the wound interface and into a secondary absorptive dressing when clinically appropriate.

This creates the possibility of a two-zone system. The first zone is the occupied wound interface, where Viniferamine Wound Hydrogel conforms intimately to wound architecture. The second is the exudate-management zone, where an appropriate secondary dressing receives and manages excess wound fluid while protecting the treatment site. These functions are complementary. The hydrogel manages the interface. The secondary dressing manages bulk fluid and external protection.

### **Potential Relevance Across Wound Types**

The occupied-interface concept could be relevant to many wounds in which irregular geometry, exudate, microbial contamination, and biofilm represent clinical challenges. These may include diabetic foot ulcers, pressure injuries, venous leg ulcers, surgical wounds, dehiscent wounds, traumatic wounds, chronic ulcers, and other complex wounds.

The principle is independent of any particular wound shape. In fact, the more irregular the wound surface becomes, the more important conformability may be. A conventional dressing can cover an irregular wound. A conformal hydrogel has the potential to occupy the irregularity itself.

### **A New Functional Property of Wound Hydrogels**

Hydrogels have traditionally been discussed in terms of moisture management, hydration, cooling, drug delivery, or patient comfort. The Viniferamine Wound Hydrogel concept suggests another possible functional property: interfacial occupation.

If demonstrated experimentally, conformity would no longer be merely a desirable handling characteristic. It could become part of the dressing system's biological function. The important measurement would therefore not be viscosity alone. We should measure wound-surface contact.

Experiments could compare Viniferamine Wound Hydrogel with conventional hydrogel or dressing interfaces using standardized irregular wound models. Imaging could quantify total contact area, void volume, fluid pooling, and the percentage of wound topography occupied by the hydrogel.

A second phase could introduce representative wound microorganisms and determine whether increasing conformity correlates with reduced initial attachment, microcolony formation, and

mature biofilm development. Such testing would allow the occupied-interface hypothesis to move from a compelling physical model to a demonstrated biological mechanism.

## Conclusion: Occupy the Space Before Biofilm Can

The most important feature of Viniferamine Wound Hydrogel may ultimately be more than its ability to conform to a wound. It may be why that conformity matters. Biofilm needs an opportunity to establish itself. Microorganisms must reach a surface, attach, multiply, and construct the protective extracellular environment that allows a mature biofilm community to develop. An irregular wound-dressing interface containing gaps, voids, exposed surfaces, and pooled exudate may provide such opportunities. Viniferamine Wound Hydrogel proposes a different architecture.

Applied in its lower-viscosity state, the hydrogel

flows into the contours of the wound. As it warms, it increases in viscosity and preserves an intimate mold of the wound surface. The result is an occupied interface designed to minimize unfilled spaces and maintain close association with the wound bed. This leads to a simple yet potentially important principle: a biofilm cannot mature without first establishing a foothold.

If an intimate wound mold can reduce the spaces and exposed interfaces that provide that foothold, physical conformity itself may become an important component of biofilm prevention. The next scientific step is therefore clear: demonstrate that intimate conformity reduces void space, then determine whether reducing that space measurably reduces microbial attachment and biofilm development. If confirmed, the bioidentical mold would represent more than a novel method of hydrogel delivery. It would represent a new strategy for engineering the wound interface itself.

## THE OCCUPIED WOUND INTERFACE

*Viniferamine® Wound Hydrogel – Exact Mold. No Gaps. Block Biofilm.*

